

PH60S-E replacement pH probe is compatible with the following pen testers:

- PH60S, PC60, PH60, PH60F

Technical Specs:


Range	0 to 14 pH	Reference	Ag/AgCl
Housing	Polycarbonate	Temp. sensor	30K Ω thermistor
Junction	Open	Electrolyte	Gel KCL
Sensor Type	LabSen® spear glass membrane	Operating temp.	32 to 122°F (0 - 50°C)

How to replace the probe?

Screw off the probe ring, unplug the probe, plug in the new replacement probe (pay attention to the probe's position), and screw on the probe ring.

Notes

1. The spear probe must be stored in the 3M KCL storage solution when not in use.
2. Never store pH probes in purified water such as RO water, tap water, distilled water, or deionized water as they will cause damage to the probe.
3. pH Probes don't last forever. Every pH probe will eventually age and fail even if you don't use it that often. The typical service life of Apera pH probes is 12-24 months depending on the frequency of usage and how well you keep it clean and properly stored. We recommend replacing your probe every 12-18 months to ensure the best accuracy.



Limited Warranty

How Long Does the Coverage Last?

Apera Instruments® (Apera) warrants the PH60S-E Replacement Probe (Product) for a period of 6 months from date of purchase by original purchaser or consumer. Proof of purchase is required for the warranty to be effective (store sales receipt for Product showing model number, payment and date of purchase). This warranty is non-transferable and terminates if the original purchaser/consumer sells or transfers the Product a third party.

What is Covered?

Apera warrants the Product against defects in material and workmanship when used in a normal manner, in accordance with Apera instruction manuals. If Apera is provided with valid proof of purchase (as defined above) and determines the Product is defective, Apera may, in its sole discretion either (a) repair the Product with new or refurbished parts, or (b) replace the Product with a new or refurbished Product.

What is NOT Covered?

This warranty does not apply to equipment, component or part that was not manufactured or sold by Apera, and shall be void if any such item is installed on a Product. Further, this warranty does not apply to replacement of items subject to normal use, wear and tear and expressly excludes:

- Cosmetic damage such as stains, scratches and dents
- Damage due to accident, improper use, negligence, careless operation or handling of Product not in accordance with Apera instruction manuals, or failure to maintain or care for Product as recommended by Apera
- Damage caused by use of parts not assembled/installed as per Apera instructions
- Damage caused by use of parts or accessories not produced or recommended by Apera
- Damage due to transportation or shipment of Product
- Product repaired or altered by parties other than Apera or its authorized agents
- Product with defaced, missing or illegible serial numbers
- Products not purchased from Apera or an Apera-authorized distributor or reseller.

How Do You Get Service?

To begin a warranty claim, have your valid proof of purchase ready and contact us at info@aperainst.de to get an RMA number.