

## AGITADOR MAGNÉTICO CON CALENTAMIENTO 380°C MSH-20D



Certificación CE y  
garantía de 2 años



Uniformidad en la  
temperatura y agitación



Incluye una barra  
magnética (30mm)



Temperatura máx.  
380 °C

### APLICACIONES

✓ SECADO DE MUESTRAS  
✓ LABORATORIO

### ESPECIFICACIONES TÉCNICAS

Spec.	Model	A. 180 × 180mm Plate (II-A) Digital "MSH-20D"
Plate Size(mm)		180× 180mm Plate-type
Material	Body	Powder Coated Aluminum
	Plate	Ceramic Coated Aluminum (Chemical / Acid Resistance, Easy to Clean)
Stirring Capa. & Motor Rating		Max. 20Lit, 9 / 4W
Temp. Uniformity		Specially Designed Heating Module : Temp. Difference 10% less
Heating Power/Consumption		600W
Stirring Bar		3cm bar-included, Usable : up to 5cm bar
Permissible Temp. & RH %		Ambient + 5 ~ 50 °C , 85% RH
Dimension & Weight		206× 307× h99mm, 3.3kg
Packing Size & Gross Weight		400× 280× h200mm, 5kg
Electrical Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz
III. Digital Package-set		(III-A) "MSH-20D" Set
Cat. No.	230V	DH.WMH03021
	120V	DH.WMH02021

+57 (1) 300 3161  
+57 316 532 6905



info@quimicompany.com.co  
comercial@quimicompany.com.co

