



## BOMBA DE VACÍO DE DIAFRAGMA

DAIHAN® 105mbar Oil Free Vacuum Pump “PU-300”, with Fine Vacuum Regulator and Moisture Trap, 20Lit/min With Prote™ Overflow Protection, Built-in Thermal Protection, Protection Class IP30.



105mbar Oil Free Vacuum Pump  
“PU-300”

- CE, CSA, and RoHS Certified
- Used for General Filtration, Microbiological Test, Suspended Solid Test, Solid-Phase Extraction, Microplate Vacuum Manifold, Air Sampling, etc.

- \* Driven by Piston without Lubricant, Regular Oil Change and Maintenance ; No Oil Pollution
- \* Prote™ Overflow Protection is Installed to Prevent a Large Amount of Liquid from being Sucked into the Pump.
- \* Equipped with Filter Cartridge in Air Inlet to Filter Particle and Moisture to Prolong the Life of Pump.

Spec.	Model	“PU-300”
Speed of Evacuation		20 Lit/min
Max. Vacuum		105 mbar
Motor Power & Max. Current		60 W, 0.3A
Motor Rotation & Noise Level		1450 rpm, 50 dB
Hose Barb		8 mm
Size & Weight		268× 135× h204 mm, 4.1kg
Electric Rqts.		220V, 50Hz
Cat. No.		DH.Pu0300

\* Other Specifications are available upon Customer's Request.

