

Microprocessor Centrifuge 1.5~50 ml Class-I Medical Device(NIDS) , Max. 4000rpm

With 6-Holes or 12& 24-Holes Fixed Angle Rotor,Digital Model, Safety Interlock, 1.5~50ml



2-Year Warranty

"Cef-D50.6" for 6 Tubes of 50ml, with 6 Holes Fixed Angle Rotor, 4000rpm



①TA50

(1) Fixed Angle Rotor ; 50ml 6Holes, "R5006"



②TA15

(2) Fixed Angle Rotor ; 10~15ml 12Holes, "R1512"



②TA15

(3) Fixed Angle Rotor ; 10~15ml 24Holes, "R1524"



(4) Swing Rotor ; 4x 50ml, 12x 1.5/2ml & 10~15ml, 16x 5~7ml "SR1550"
Included ⑤,⑥,⑦ / Option ④



①TA50 ②TA15

Optional Tube Adapter

- **DAIHAN Designed** ■ **Durable and Maintenance Free - Brushless DC Motor**
- **Ideal for Biology, Medicine, Pharmacology, Laboratory, etc. Various Application**
- **with 6-Holes or 12- & 24-Holes Fixed Angle Rotor for 1.5ml ~ 50ml Tubes ; On Adapters Adjustable 1.5/2ml, 5~7ml, 10~15ml, and 50ml**
- **Automatic Quick Self-Cooling by Fan** ■ **Interchangeable Rotor**
- * Compact Design with Automatic Door Release. * Easy to Use Controls and Display
- * Self-diagnostic Function : Lid Opening, Unbalance, Over-Heat/ Current/ Voltage/ Speed
- * Auto-balancing System for Extend Motor Life and Reduce Noise
- * Program Memory Function : 9Set. * Precision Digital Timer, Set to 10Sec. / Max 30 Min
- * Pushing the Function Key for Short Runs * Rapid Deceleration due to Automatic Brake System

Specifications

| Spec. | | Model | "Cef-D50.6" |
|-------------------------------|----------|---|--|
| | | | 6 Holes - 50ml or 12/24 Holes - 1.5~50ml |
| Fixed Angle Rotor&Capa. | Included | | (1) 6Hole ; 50ml Tubes |
| | Option | | (2) 12Hole ; (1.5~15)ml Tubes, (3) 24Hole ; (1.5~15)ml Tubes |
| Swing Rotor & Capa. | Option | | (4) SR1550 ; ⑤ 4x 1Hole - 50ml, ⑥ 4x 3Hole - 10~15ml, ⑦ 4x 4Hole : 5~7ml, ⑧+⑨ 4x 3Hole : 1.5/2.0ml |
| | | | (1) 6 Tubes ; 50ml |
| Useful Tube size & Acceptance | | When Option Fixed Angle Rotor | When Option Swing Rotor ; (4) SR1550 |
| | | (2) 12 Tubes ; 1.5~15ml (3) 24 Tubes ; 1.5~15ml | ⑤ 4 Tubes : 50ml, ⑥ 12 Tubes : 10~15ml ⑦ 16 Tubes : 5~7ml ⑧+⑨ 12 Tubes : 1.5/2.0ml |
| Adapter and Size | Included | | ① 6eax 50ml Tubes, Size ; id30.7x L104mm |
| | Option | When Option Fixed Angle Rotor ② 12/24eax 10~15ml Tubes Size ; id17.5x L95mm ③ 12/24eax 5~7ml Tubes Size ; id13.5x L69mm ④ 12/24eax 1.5/2ml Tubes Size ; id11x L42mm | When Option Swing Rotor ; ⑤ 4x 1Hole - 50ml Tubes ⑥ 4x 3Hole - 10~15ml Tubes ⑦ 4x 4Hole - 5~7ml Tubes ⑧+⑨ 4x 3Hole - 1.5/2.0ml Tubes |
| Speed Range (rpm) & RCF | | 600~4000rpm & 2,325x g (Swing Rotor ; 2,683x g) | |
| Motor | | Brushless DC Motor | |
| Timer | | 0~30min | |
| Dimension & Weight | | 390x 505x 290mm, 19.8kg, Gross. 21.8kg | |
| Electric Reqts. | | AC 120V, 60Hz | |

Spare Rotor and Tube Adapter

- DH.CefR5006** (1) 50ml 6Holes Angle Rotor "R5006" with 6eax 50ml Tube Adapter PP "TA50", idΦ 30.7x L104mm, for "Cef-50.6" & "Cef-D50.6"
- DH.CefA050** ① 50ml Tube Adapter PP "TA50", idΦ 30.7x L104mm, Fit to Rotor "R5006", for Centrifuge "Cef-50.6" & "Cef-D50.6"
- DH.CefA015** ② 10~15ml Tube Adapter PP "TA15", idΦ 17.5x L98mm, for Centrifuge "Cef-6" & "Cef-8", and Rotor "R1512" & "R1524"

Optional Rotor and Tube Adapter

- DH.CefR1512** (2) 10~15ml 12Holes Angle Rotor "R1512", with 12eax 10/15ml Tube Adapter "TA15", idΦ 17.5x L98mm, for Centrifuge "Cef-50.6" & "Cef-D50.6"
- DH.CefR1524** (3) 10~15ml 24Holes Angle Rotor "R1524", with 24eax 10/15ml Tube Adapter "TA15", idΦ 17.5x L98mm, for Centrifuge "Cef-50.6" & "Cef-D50.6"
- DH.CefSR1550** (4) Swing Rotor "SR1550", with Tube Adapter PP 1Holex 4eax 50ml "TA5001", 3Holex 4eax 1.5/2ml & 10~15ml "TA1503" and 4Holex 4eax 5~7ml "TA704", for Centrifuge "Cef-D50.6"
- DH.CefA007** ③ 5~7ml Tube Adapter PP "TA7", idΦ 13.5x L69mm, Fit to "TA15", Tube Adapter, for Centrifuge "Cef-6", "Cef-8", "Cef-50.6", and "Cef-D50.6"
- DH.CefA002** ④ 1.5/2ml Tube Adapter PP "TA2", idΦ 11x L42mm, Fit to "TA15", "TA1503" Tube Adapter, for Centrifuge "Cef-6", "Cef-8", "Cef-50.6", and "Cef-D50.6"
- DH.CefA5001** ⑤ 50ml Tube Adapter PP "TA5001", Fit to Rotor "SR1550", for Centrifuge "Cef-D50.6"
- DH.CefA1503** ⑥ 10~15ml 3Hole Tube Adapter PP "TA1503", Fit to Rotor "SR1550", for Centrifuge "Cef-D50.6"
- DH.CefA704** ⑦ 5~7ml 4Hole Tube Adapter PP "TA704", Fit to Rotor "SR1550", for Centrifuge "Cef-D50.6"

+57 (1) 300 3161
+57 316 532 6905

info@quimicompany.com.co
comercial@quimicompany.com.co

