




## MICROSCOPIO BINOCULAR SERIE BM33B



Binocular



Cabeza giratoria  
360°



Brillo ajustable




Iluminación  
Luz fría



Campo  
brillante



Alta  
resolución



### ESPECIFICACIONES TÉCNICAS

Spec.	Model	BM33B
Optical System		Infinite Optical System
Viewing head		Seidentopf Binocular Viewing Head Inclined at 30°, interpupillary 55-75mm, 360° Rotation
Nosepiece		Backward Quadruple revolvingnosepiece
Eyepiece		Wide Field Eyepiece WF10x/18 mm
Objectives		Achromatic objective: 4x/0.10 10x/0.25 40x/0.65, spring loaded 100x/1.25, spring loaded, oil
Stage		Double Layer Mechanical Stage Size 125x125mm/60x40mm
Focusing		Coaxial Coarse and Fine Adjustment, Fine Division 0.002mm, Moving Range 14mm
Condenser		NA 1.25 Abbe Condenser with Iris Diaphragm
Illumination		High brightness energy-saving LED illumination for long life, brightness adjustable
Supplied with		Dust cover, spare fuse, blue, green and yellow filters
Supplied with		180 x 352 x 220 mm, net weight: 4 kg
Power Consumption		110-220V, 50/60Hz

- ✓ Universidades
- ✓ Laboratorios
- ✓ Veterinarias
- ✓ Médicas