



MICROSCOPIO TRINOCULAR SERIE IC2100T



Trinocular



Cabeza giratoria
360°



Iluminación Kohler
opcional



Configuración contraste de fase,
campo oscuro y polarización básica



ESPECIFICACIONES TÉCNICAS

Spec.	Model	IC2100T
Optical System		Infinite Optical System
Viewing head		Seidentopf Trinocular Viewing Head Inclined at 30°, interpupillary 47-75mm, 360° Rotation
Nosepiece		Backward Quadruple Nosepiece
Eyepiece		Wide Field Eyepiece EW10x/20
Objectives		Infinite semi-plan achromatic objective: 4x/0.10 10x/0.25 40x/0.65, spring loaded 100x/1.25, spring loaded, oil
Stage		Double Layer Mechanical Stage 150X140mm, Moving Range 75X50mm
Focusing		Coaxial Coarse and Fine Adjustment, Fine Division 0.002mm, Coarse Stroke 37.7mm per Rotation, Fine Stroke 0.2mm per Rotation, Moving Range 25mm
Condenser		NA 1.25 Sliding-in Centerable Condenser
Illumination		LED illumination for long life, brightness adjustable
Supplied with		Dust cover, spare fuse, blue & green filters
Power Consumption		100-240V, 50/60Hz

- ✓ Universidades
- ✓ Laboratorios
- ✓ Industria
- ✓ Ciencias biomédicas