

## CARACTERÍSTICAS



Dispositivo médico tipo I



Equilibrio automático



Velocidad estable



Amplia gama de accesorios

## ESPECIFICACIONES TÉCNICAS

Spec.		Model	I. 15ml "Cef-6" 6 Holes - 1.5~15ml	II. 15ml "Cef-8" 8 Holes - 1.5~15ml	III. 50ml "Cef-50.6" 6 Holes - 50ml or 12/24 Holes - 1.5~15ml
Fixed Angle Rotor&Capa.	Included		6Hole ; (1.5~15)ml Tubes	8Hole ; (1.5~15)ml Tubes	(1) 6Hole ; 50ml Tubes
	Option				(2) 12Hole ; (1.5~15)ml Tubes (3) 24Hole ; (1.5~15)ml Tubes
Useful Tube size & Acceptance			6 Tubes ; 10~15ml, 1.5~7ml on Adapters	8 Tubes ; 10~15ml, 1.5~7ml on Adapters	(1) 6 Tubes ; 50ml When option Rotor ; (2) 12 Tubes ; 10~15ml (3) 24 Tubes ; 10~15ml 1.5~7ml on Adapters
Adapter and Size	Included		② 6eax 10~15ml Tubes Size ; id17.5x L95mm	② 8eax 10~15ml Tubes Size ; id17.5x L95mm	① 6eax 50ml Tubes Size ; id30.7x L104mm
	Option		③ 6eax 5~7ml Tubes Size ; id13.5x L69mm ④ 6eax 1.5/2ml Tubes Size ; id11x L42mm	③ 8eax 5~7ml Tubes Size ; id13.5x L69mm ④ 8eax 1.5/2ml Tubes Size ; id11x L42mm	② 12/24eax 10~15ml Tubes Size ; id17.5x L95mm ③ 12/24eax 5~7ml Tubes Size ; id13.5x L69mm ④ 12/24eax 1.5/2ml Tubes Size ; id11x L42mm
Speed Range (rpm)			300~3000rpm at 60Hz (300~2500rpm at 50Hz)	300~6000rpm	600~4000rpm
R.C.F			1,066x g	4,265x g	2,325x g
Display			None	Tachometer for Speed & R.C.F	None
Timer			0~30min Analog Timer		
Motor			Induction Motor	DC Motor	
Safety Manual Brake			None	Yes	None
Dimension & Weight			280x 280x h255mm, Net. 3.5kg, Gross. 4kg	280x 280x h265mm, Net. 4.5kg, Gross. 5kg	300x 300x h315mm, Net. 8.5kg, Gross. 9kg
Electric Reqts			AC 230V, 50/60Hz or AC 120V, 60Hz		
Cat. No.	230V		DH.Cef0106	DH.Cef0108	DH.Cef0110

+57 (1) 300 3161  
+57 316 532 6905



info@quimicompany.com.co  
comercial@quimicompany.com.co

