



QUIMICOMPANY

EQUIPOS E INSUMOS PARA LABORATORIO

APERA
INSTRUMENTS

PH PEN TESTER; SUITABLE FOR TEST TUBE PH TESTING

REF: SX610



- Easy SINGLE-Button Operation, Smallest Pen-Style design with a replaceable probe (good for test tube measurements)
- AUTO Calibration with auto buffer recognition; Auto Temperature Compensation (ATC) ensures accuracy in variant temperatures 0-60°C (32-140°F)
- COMPLETE KIT with READY-TO-USE calibration solution, KCl storage solution, and batteries in a portable carrying case
- Intelligent functions: automatic power-off, low battery alarm, and pH probe invalidation reminder
- AUTO RECOGNITION and REMINDER of stable values improves consistency, saving you the trouble of guessing when to manually hold measurement

Is featured with reliable measuring performance with a single-button operation in a slim body. A cost-effective choice for your regular pH testing, including hydroponics, water treatment, education, pools& spas, and etc.

Kit contents



The SX610 pH Pen Tester/ Meter comes with a complete kit of premixed calibration solutions (7.00 pH), storage solutions (3M KCL), CR2032 lithium batteries all in a portable carrying case. IP57 waterproof and dustproof floats on water, so you don't have to worry if it falls into the water by accident.

☎ (601) 300 3161 + (57) 316 532 6905

🌐 www.quimicompany.com.co

✉ info@quimicompany.com.co