

AGITADOR SHAKER MaXshake™ OB2



LCD display



Alarmas audibles cuando ocurre un error



Control de velocidad ajustable



Control de velocidad por microprocesador

ESPECIFICACIONES TÉCNICAS

Spec.	Model	A. "MaXshake™ OB2" Orbital Motion
Shaking Platform Size		255 × 255 mm
Shaking Motion		Orbital Motion
Speed Range & Max. Load		10 ~ 300 rpm, 10kg
Amplitude of Rotation (mm)		10 mm
Motor		Brushless DC Motor
Permissible Temp. & Moisture		Ambient +5 ~ 50 °C 80%
Timer & Alarm		99hr 59min (with Continuous run.), Error status & Timer-end
Display		Large Size Digital LCD with Back-Light Function
Controller		Digital Fuzzy Controller with Touch Button, Jog-dial
Program Function		10 Steps, 6 Programs
Platform		Universal Platform for various Shapes & Size of Flask
Dimension & Weight		376× 375× h138 mm(Base Platform : 255× 255 mm), 13kg
Packing Dimension & Weight		345× 420× h205 mm, 13.5kg
Safety Circuit		Overload Protector for Motor
Power Consumption		45 W
Electrical Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz

Cat. No. 120V

DH.WSO00301

+57 (1) 300 3161
+57 316 532 6905



info@quimicompany.com.co
comercial@quimicompany.com.co

