


MULTI HEATING MANTLE WITH 6 PLACES 300ML

REF:DH.WHM12294

DAIHAN® Aluminium-case Flask Heating Mantle, (1) Basic & (2) Stirring-type, 450 °C, 50~20,000ml
With Built-in Temp Controller, with/without Mag-stir Speed Control, with Certi. & Traceability



- Ideal for DURAN / Kimax / Pyrex & Iwaki-brand Round Flasks
- With Certificate & Traceability : Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Ideal for Repetitive Extraction (Kjeldahl, Soxhlet and other Extraction Procedure), Reflux and Distillation Laboratory Applications
- Analog Controller Built-in System
- Rods and Clamps for Rods Connection
- Highly-processed Aluminum Case
- Nickel Chrome Heating Element - Coil Type
- Without Clamp & Clamp Holder for Glassware ; Order Separately

Premium Analog Multi Heating Mantle

Highly-processed Aluminum Case
Nickel Chrome Heating Element

Spare Support Rod & Clamp Set

DH.WHM60250 Support Rod & Clamp Set, "MHR03", for 6Places 250-/300-ml, 2×750-/1×900-mm Rod and 2×Clamp for Rod

Remarks : DU. : DURAN Ki. : Kimax Iw. : Iwaki Py. : Pyrex

Technical specifications

Parameters		Description
Electrical requirement		120V/60Hz
Flask	Capa. (ml)	250/300
	Dia. Max. Φ (mm)	83 ~ 86
	Bottom Style	Round
Watts(W)		840
Inside Depth (mm)		50
Dimension w×d×h (mm)		800×240×110

☎ (601) 300 3161 +(57) 316 532 6905

🌐 www.quimicompany.com.co

✉ info@quimicompany.com.co