

HEATING INCUBATOR - 42 L

REF: BPH-9042

LCD Programmable controller (with timing function)

Features

- Polished stainless-steel chamber, semicircular arcs at corners for easy cleaning, and the space between the shelves in the chamber is adjustable.
- Auto-controller of fan speed to prevent damage to the samples.
- Large LCD screen to display more data at same time.
- Self-check function easy to identify problems.
- Over temperature and temperature difference alarms.
- There is a 25mm instruction connection hole on the left side of the chamber for easy testing operation and temperature measurement.
- Programmable controller: 7 periods 63 steps, 0 ~ 5999mins for each periods, fan speed 0 to 100% adjustable.
- Independent audible and visible temperature-limiting alarm system ensures experiments run safely.
- RS485 connector can connect computer and printer to record the parameters and the variations of temperature.(Option)



Technical specifications

Parameters	Description
Electrical Requirement	110V 60Hz
Controller	LCD Programmable controller
Temperature Range	RT+5°C ~80°C
Display Resolution	0.1°C
Temperature Stability	±0.2°C
Temperature Uniformity	±1.5(at 37°C)
Ambient Temperature	+5~35°C
Power Consumption	250W
Volume	42L
Interior Dimension(W×D×H)mm	350x300x400
Exterior Dimension(W×D×H)mm	505x635x600
Shelves	2(pcs)
Timing Range	0~5999min
Remark	Model BPH-9042 without observation window

Options



USB data collect



RS485 connector



Printer



UV Sterilizer

(601) 300 3161 +(57) 316 532 6905

www.quimicompany.com.co

info@quimicompany.com.co