

BASIC STEREO MICROSCOPE - 90° ROTATION

REF: ST57-3C



Descriptions

The ST57 series stereo microscope is an extremely robust, stable stereo microscope which is easy to use, it is ideal for all conventional applications in schools, workshops and training companies.

The whole range is assembled according to the best technical standards. The optical system, thanks to its transparency, quality and depth, is at the top of its category. The series is composed of 15 models, each one designed to satisfy multiple needs.

Objective Specifications

Eyepiece	Magnification	1x	2x	3x	4x
WF10x	Total magnification	10x	20x	30x	40x
	Field of view mm	∅ 20	∅ 10	∅ 6.7	∅ 5
Working Distance		57cm	57cm	57cm	57cm

BASIC STEREO MICROSCOPE - 90° ROTATION

REF: ST57-3C



Technical specifications

Parameters	Description
Interpupillary distance	55~75
One tube with diopter adjustment	+/-5 Dio
Illumination	Incident/Transmitted Light
Optional accessories	Inserted Objective 1X, 3X, 4X, 6X
Head and Body	Inclined Viewing Un-rotatable Body
Objective	2X & 4X 90° Rotated
W.D.	57mm
Eyepiece	WF10X/20mm
Electrical requirement	110V or 220V/12V 10W