

INTELLIGENT VISCOMETER 20-100.000mPa.S

REF: NTV-E1



Features

- Metal lifter
- Build-in temperature probe
- New designed durable shaft
- Auto switch between dynamic viscosity and kinematic viscosity
- Variety of viscosity units
- Continuous viscosity value displayed
- Over range alarm
- Linear calibration by computers
- Power supply: 1 00V-240V, good anti-interference

Optional accessories

- Thermostatic water bath (for viscosity measurement)
- Ultra-low viscosity adapter
- Standard oil
- Micro printer
- External circulatory glass beaker

Applications

For low and medium viscosity measurement

Technical specifications

Parameters	Description
Range (mPa.s)	20-100.000
Rotate Speed (RPM)	6,12,30,60
Spindles	LV 1, 2, 3, 4 Spindle (LV O is optional)
Accuracy	± 1.0%(FSR)
Repeatability	± 0.2%(FSR)