



QUIMICOMPANY

EQUIPOS E INSUMOS PARA LABORATORIO

Denmark
SCIENTIFIC

INTELLIGENT VISCOMETER 20-2.000.000mPa.S

REF: NDJ-8ST



Features

- Anti-static shell and PC lifter
- RTD temp probe
- Auto switch between dynamic viscosity and kinematic viscosity
- Variety of viscosity units
- Accompanied with operation manual and viscosity calibration
- Continuous viscosity value displayed
- Power supply: 1 00V-240V, good antiinterference

Optional accessories

- Thermostatic water bath (for viscosity measurement)
- Ultra-low viscosity adapter
- Standard oil
- Micro printer
- External circulatory glass beaker

Applications

For low and medium viscosity measurement

Technical specifications

Parameters	Description
Range (mPa.s)	20-2.000.000
Rotate Speed (RPM)	0.3,0.6,1.5, 3, 6, 12, 30, 60
Spindles	LVI, 2, 3, 4 Spindle (LV O is optional)
Accuracy	± 1.0%(FSR)
Repeatability	± 0.5%(FSR)

(601) 300 3161 +(57) 316 532 6905

www.quimicompany.com.co

info@quimicompany.com.co