

SIEVER WITH TOUCH SCREEN REMOTE CONTROLLER

REF: BNVSS-C200




Features

- Particle analysis
- Widely used in the field of pharmaceutical, metallurgy, food, cosmetics
- Advanced integrated circuit control
- Easy to operate
- Electronic control of amplitude and frequency
- The sieves are chosen by users

Technical specifications

Parameters	Description
Operation range	≤325mesh
Suitable sieve diameter	Φ200mm
Sieve height	50mm
Max Layers	1-6 layers
Vibration frequency	3000 times/min, 6000 times/ min
Amplitude selection	0 - 3mm continuous adjustment
Shaking modes	a).Subtle shaking
	b).Intermittent shaking, amplitude can be adjusted within the range of 1 -2 and 2 -3, increment: 0.2±0.1mm
	c).Continuous shaking, amplitude can be adjusted within the range of 1 -2 and 2 -3, increment: 0.2±0.1mm
	d).Three modes can be precisely timed
	e).Three modes can be memorized even when the equipment is shut off
Electrical supply	AC 110V 50HZ
Dimensions	Shaking unit: Φ 300×210mm
	Control unit: 280×200×140mm
Weight	25kg



SIEVER WITH TOUCH SCREEN REMOTE CONTROLLER

REF: BNVSS-C200

Test Sieves

Mesh Size	Wire Dia(mm)	Aperture Size(mm)
4	1	5
5	1	4
6	0.8	3.35
7	0.8	3
8	0.8	2.36
10	0.63	2
12	0.4	1.7
14	0.4	1.4
16	0.4	1.18
18	0.315	1.12
20	0.315	1
24	0.28	0.8
26	0.28	0.71
28	0.28	0.63
30	0.25	0.6
32	0.224	0.56
35	0.224	0.5
40	0.18	0.45
45	0.14	0.425
50	0.14	0.355
55	0.14	0.315
60	0.14	0.28
65	0.14	0.25
70	0.125	0.224
75	0.125	0.212
80	0.125	0.2
85	0.125	0.18
90	0.1	0.16

Mesh Size	Wire Dia(mm)	Aperture Size(mm)
100	0.1	0.15
120	0.09	0.125
130	0.071	0.12
140	0.071	0.112
150	0.063	0.106
160	0.056	0.1
180	0.063	0.08
190	0.056	0.078
200	0.05	0.075
220	0.045	0.071
240	0.05	0.06
250	0.045	0.056
260	0.04	0.053
280	0.04	0.05
300	0.036	0.048
320	0.036	0.047
325	0.036	0.045
350	0.032	0.04
360	0.032	0.038
400	0.028	0.036
500	0.025	0.028
600	0.02	0.02
800		0.018
900		0.015
1000		0.013
1100		0.012
1600		0.008
2000		0.0065



- Material: Stainless steel
- Frame diameter: 200mm / 300mm
- Height: 50mm / 60mm
- Standards: GB/T5330-2003
- Available: Wire mesh sieves and Perforated plate sieves
- Sieve accessories: Pans, lids